

QUALITY Care is Only Fair



**MEHTA
SURGICARE
PVT. LTD.**

ABOUT US

Mehta Surgicare Pvt. Ltd. is a well-established company since 25 YEARS. Company is expertise in manufacturing surgical dressing products under revised schedule 'M as well as under schedule' T of Drug and Cosmetics act 1945.


The company is having well equipped with latest technology, machineries established to manufacture Surgical Adhesive Tape, Belladonna Perforated Plaster, Capsicum–Rapid Pain Reliever, Capsicum Belladonna Perforated Plaster, Cornex- corn dressing, Boilex-Boil dressing, Medicated First Aid Dressing, (Non-Woven Woven, Cloth as well as wash proof) Bandages (Cotton Crap Bandage & Elastic Crap Bandage of various size), Medifix-Cannula Fixator, Micropourus Non-woven Surgical Paper Tape, 3 Ply, NK95 etc. Masks, Medical Devices & Hospital Division disposables etc.

To ensure Company continue to deliver on our commitments, Company have developed High-Tech. Laboratory governed by skillful technical staff which processes the quality products.

The foundation of guiding value are:

* Quality, Fair & Transparent Business Association, Customer Satisfaction, Reliability Consistency, Innovation. Mehta Surgicare Pvt. Ltd. remain steadfast to its goal of providing better quality products to valuable customer. Mehta Surgicare Pvt. Ltd. also has adopted an international quality management system and have the ISO 9001:2015, ISO 13485:2003 certifications along with self-product certification awarded by CE mark.

Mehta Surgicare Pvt. Ltd. is awarded with GMP (Good Manufacturing Practice) certificate from Govt. of Gujarat (India). Company.



**MEHTA
SURGICARE
PVT. LTD.**

Q U I L I T Y Care is Only Fair




mehta's
Medistrip[®]
MEDICATED FIRST AID DRESSING

Relieves: Mehta's Medicated strip first aid dressing pad gives antiseptic protection for children, players, sportsman, athlete shoe bite, house wife working with kitchenware like knife, pillar, shredder and seizer.

Active Constituents:
Each pad contains Haldi Powder (Curcuma longa) 0.12% w/w

Application / Direction: Remove the plastic backing strip and the position the felt turmeric coating over the wound so that medicated centre cover the wound or injured area. Apply as & when require. Easy to apply and remove. Perforation for ventilation. Contains proven antiseptic. Apply to clean and dry skin. Hygienically packed. Apply to clean and dry skin





**MEHTA
SURGICARE
PVT. LTD.**

QUALITY Care is Only Fair



BELLADONNA PERFORATED PLASTER

- SCIATICA
- LUMBAGO
- RHEUMATISM
- DISTORTIONS DUE TO WHIPLASH
- BOILS
- MUSCLE STIFFNESS
- SWELLING PAIN
- MUMPS
- MUSCULAR PAINS
- BACKACHE

SIZE :

30 Box / Case (30 Inner x 100) 10 Pouch x 10 Pcs.

Application / Direction : Belladonna Perforated Plaster to be applied externally. Clean the affected area(s) with soap and water. Dry the skin thoroughly. Remove the protective covering from the face of the plaster and apply to affected clean and dry skin by pressing plaster firmly on to the skin. The application can stay on skin for 2 days and then removed. If area is still painful, then apply new plaster. Belladonna Perforated Plaster can be used until relief from pain, use as often as needed.

CAPSICUM + BELLADONNA

PERFORATED PLASTER


- MUMPS
- BACKACHE
- SWELLING PAIN
- DISTORTIONS DUE TO WHIPLASH
- BOILS
- MUSCULAR PAINS
- MUSCLE STIFFNESS
- SCIATICA
- RHEUMATISM
- LUMBAGO

SIZE :

30 Box / Case (30 Inner x 100) 10 Pouch x 10 Pcs.



Application / Direction: CAPSICUM + BELADONNA PERFORATED PLASTER to be applied externally. Clean the affected area(s) with soap and water. Dry the skin thoroughly. Remove the protective covering from the face of the plaster and apply to affected clean and dry skin by pressing plaster firmly on to the skin. The application can stay on skin for 2 days and then removed. If area is still painful, then apply new plaster. CAPSICUM + BELADONNA PERFORATED PLASTER can be used until relief from pain, use as often as needed.



**MEHTA
SURGICARE
PVT. LTD.**

QUALITY Care is Only Fair



mehta's
CAPSICUM
RAPID PAIN RELIEVER

SCIATICA **LUMBAGO** **RHEUMATISM**

DISTORTIONS DUE TO WHIPLASH

MUSCLE STIFFNESS

SIZE :

30 Box / Case (30 Inner x 100) 10 Pouch x 10 Pcs.

Application / Direction: MEHTA'S CAPSICUM RAPID PAIN RELIEVER to be applied externally. Clean the affected area(s) with soap and water. Dry the skin thoroughly. Remove the protective covering from the face of the plaster and apply to affected clean and dry skin by pressing plaster firmly on to the skin. The application can stay on skin for 2 days and then removed. If area is still painful, then apply new plaster. MEHTA'S CAPSICUM RAPID PAIN RELIEVER - PERFORATED ADHESIVE PLASTER can be used until relief from pain, use as often as needed.

CORNEXTM
Corn Caps


MEHTA'S CORN-X CAPS ACTS TO RELIEVE AND REMOVE HARD CORNS AS WELL AS WARTS.

SIZE :

27 Box / 50 Pcs. = 1350 Pcs. / Case | 4 Pcs. Per Card



Application / Direction: Ensure feet are clean and dry. Remove plastic backing strip of Corn-x Cap and position the felt ring over the corn/wart. Firmly fix position with adhesive strips, when properly applied, it should not slip or move. Change the Corn-X caps every other day until necessary. A maximum of four Corn-x should be used on each application.




**MEHTA
SURGICARE
PVT. LTD.**

QUALITY Care is Only Fair

BOILEX™

MEDICATED BOIL PLASTER

**MEHTA'S BOILEX - BOIL DRESSING MEDICATED
PLASTER RELIEVE AND REMOVE BOILS.**



Application / Direction: Ensure feet are clean and dry. Remove the plastic backing strip and position the felt ring over the boil, so that the medicated centre cover the boil. When applied must NOT move or slip. Remove Boilex dressing every twelve (12) hours, clean affected area and replace with a new dressing. If boil is not cured within two days, repeat process using Boilex dressing.

WASH PROOF

Medistrip

MEDICATED FIRST AID DRESSING



**FIRST AID DRESSING FOR
CUTS AND WOUNDS**


Relieves: Mehta's Medicated strip first aid dressing pad gives antiseptic protection for children, players, sportsman, athlete shoe bite, house wife working with kitchenware like knife, pillar, shredder and seizer.

Active Constituents:

Each pad contains Haldi Powder (Curcuma longa) 0.12% w/w



Application / Direction: Remove the plastic backing strip and the position the felt turmeric coating over the wound so that medicated centre cover the wound or injured area. Apply as & when require. Easy to apply and remove. Perforation for ventilation. Contains proven antiseptic. Apply to clean and dry skin. Hygienically packed. Apply to clean and dry skin



**MEHTA
SURGICARE
PVT. LTD.**

QUALITY Care is Only Fair

mehta's
Radiplast
MEDICATED
FIRST AID DRESSING

**FIRST AID DRESSING FOR
CUTS AND WOUNDS**



Indication: RADIPLAST (Ayurvedic Medicine), protect the wounds against the infections, Heal the wounds, it prevent the inflammation of the wounds. It contain proven antiseptic, easy to apply and remove.

Application / Direction: Remove the plastic backing and the position the felt turmeric coating over the wound so that medicated centre cover the wound or injured area. Apply as & when require. Easy to apply and remove. Perforated for ventilation. Contains proven antiseptic, Hygienically packed. Apply to clean and dry skin.

mehta's

ADHESIVE TAPE


ACTIVE CONSTITUENTS: MADE OF STRONG AND FLEXIBLE FIBRES THAT CONTRACT AND EXPAND WITH MOVEMENT. MEDICAL CLOTH TAPES ARE USUALLY STRONG, ELASTIC IN NATURE AND ABLE TO HOLD HEAVIER TUBING AND THICK DRESSINGS.

SIZE :

1.25 cm x 1 Mtr. (1350 Box / Case)	30 Box x 45 Pcs.
2.25 cm x 1 Mtr. (1440 Box / Case)	48 Box x 30 Pcs.
2.25 cm x 5 Mtr. (420 Box / Case)	21 Box x 20 Pcs.
5 cm x 5 Mtr. (210 Box / Case)	21 Box x 10 Pcs.
7.5 cm x 5 Mtr. (420 Box / Case)	15 Box x 30 Pcs.



Application / Direction: Apply to clean and dry place. Mehta's Adhesive tapes are used in first aid to secure the bandage or a wound dressing in their place of application. Mehta's Adhesive tapes are also used to hold IV lines, oxygen tubing or splints in position. Waterproof Surgical tapes may be used to hold together a gaping wound while the wound is prepared for surgical sutures or it may be used in place of stitches to keep the wound closed while it heals.



**MEHTA
SURGICARE
PVT. LTD.**

QUALITY Care is Only Fair

MEHTAPORE
MICROPOROUS SURGICAL TAPE



ACTIVE CONSTITUENTS : MICROPOROUS SURGICAL PAPER TAPE IS MADE UP OF NON-WOVEN FABRIC, EASILY REMOVABLE AND HIGH TEARING STRENGTH. UNIFORMLY COATED WITH HYPOALLERGENIC, NON-TOXIC, NON-IRRITANT ADHESIVE

SIZE :

1.25cm x 9 Mtr.	(60 Box / Case)	24 Pcs. Per Box
2.5 cm x 9 Mtr.	(60 Box / Case)	12 Pcs. Per Box
5 cm x 9 Mtr.	(60 Box / Case)	06 Pcs. Per Box
7.5 cm x 9 Mtr.	(60 Box / Case)	04 Pcs. Per Box
1.25 cm x 5 Mtr.	(60 Box / Case)	24 Pcs. Per Box
2.5 cm x 5 Mtr.	(60 Box / Case)	12 Pcs. Per Box
5 cm x 5 Mtr.	(60 Box / Case)	06 Pcs. Per Box
7.5 cm x 5 Mtr.	(60 Box / Case)	04 Pcs. Per Box

Application / Direction: Used for fixation of dressing, tubing, post-surgical care and pathology lab etc. on patients. Uniformly coated with hypoallergenic, non-toxic, non-irritant adhesive. Easy to peel off without causing discomfort or skin irritation. High porosity allows skin to breathe, and high resistance to water and vapour prevents the tape from peeling off. Suitable on sensitive body parts, for infants and adults, Highly Breathable and Water Resistant.

mehta's

COTTON CREPE BANDAGE WITH FAST EDGES


ACTIVE CONSTITUENTS : IT CONSISTS OF A WOVEN FABRIC, FLEXIBLE IN WRAP, IN ONE CONTINUOUS LENGTH WITHOUT JOINTS WITH FAST EDGES. THIS COMES ALONG WITH CLIPS TO TIE THE BANDAGE. THE WRAP THREADS CONSIST OF TWO-FOLD COTTON THREADS CREPE TWISTED AND WELT THREADS THAT ARE MADE OF COTTON. THE FABRIC IS CLEAN AND FREE FROM WEAVING DEFECTS AND OTHER IMPURITIES.

SIZE :

8 CM x 4 Mtr. (120 Contr.) 1 Roll Per Contr.
10 CM x 4 Mtr. (120 Contr.) 1 Roll Per Contr.



Application / Direction: Mehta's Cotton Crepe Bandage used for various veins, strains, sprains, painful joints and as a pressure bandage in burns and skin grafts. Wash carefully in hot soapy water. Rinse in fresh tap water. Squeeze gently and allow to dry without stretching or wringing. On no account should it be ironed.



**MEHTA
SURGICARE
PVT. LTD.**

QUALITY Care is Only Fair

mehta's
**ELASTIC ADHESIVE
BANDAGE**

WITH FAST EDGES



SIZE :

8 cm x 4 Mtr. 1 Roll Per Contr
10 cm x 4 Mtr. 1 Roll Per Contr


Application / Direction: *Elastic Adhesive Bandage with fast edges* - flash coloured elastic with porous adhesive bandage has good aesthetic appeal due to woven fast edges, thick Fabric, allows Free muscle & Joint Movement, allow skin breathing, acts as a mild antiseptic, astringent & anti-inflammatory, long protection & long storage life.

ELASTIC ADHESIVE BANDAGE B.P.

Medifix
CANNULA FIXATOR



Application / Direction: ELASTIC ADHESIVE BANDAGE B.P. MEDIFIX CANNULA FIXATOR are Used for fixation of dressing, -surgical care and pathology lab etc. on patients. Uniformly coated with hypoallergenic, non-toxic, non-irritant adhesive. Easy to peel off without causing discomfort. High porosity allows skin to breathe, prevents the tape from peeling off. Suitable on sensitive body parts, for infants and adults, Highly Breathable.



**MEHTA
SURGICARE
PVT. LTD.**

MEHTA SURGICARE PVT. LTD.

G-523, Lodhika, G.I.D.C. Metoda-360021 Dist. Rajkot (Guj.) India.

E-mail : info@mehtaherbals.com • Web : www.mehtaherbals.com